


Cath-Secure® Lock Overview

- Low profile
- No plastic parts, resulting in lower cost to end user
- Long securing tab with adhesive back
- Able to accommodate most sizes of tubing
- Adhesive tab to prevent movement of tube
- Velcro base to prevent kinking of tube
- Easy to remove with or without alcohol
- DEHP Free, Hypoallergenic, Not Made With Natural Rubber Latex



Item No.	Description	Qty.	HCPCS	Measurements
5445-1	Foley and multi-purpose tube securement (Velcro anchor)	40/ctn - 8 ctn/cs - 320 ea/cs	A4333	3" tab & 3.12" x 1.8" base

Statlock® Lock Overview

- Single application product, higher cost to facilities
- Moving parts, plastic clip may catch on clothing/linens
- Only intended for use with Foley catheters
- Skin protectant prep pad used for application
- Alcohol needed for removal
- Foley must be inserted before placed on skin

* Plastic clip has the potential to separate from base
 Swivel base may cause tube to kink, restricting flow
 Plastic clip may not close properly

*MAUDE Reports: January 1st, 2018 - July 31st, 2019
<https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfMAUDE/search.CFM>



Manufacturer	Product	Item No.	UOM	Size	Components
M.C. Johnson Co., Inc.	Cath-Secure® Lock	5445-1	40/box	Adult	Made in the USA
Bard Medical	Statlock® Foley Stabilization Device	FOL0102	25/box	Adult	Assembled in Mexico