

# ALGINATE DRESSING

## Features and Benefits

- Non-woven
- Non-adhesive
- Highly absorbent
- Needled for extra comfort
- Removes in one piece after gelling
- Soft and comfortable



## Description

The sterile, absorbent calcium alginate dressing for use on moderate to heavily draining wounds. It may be appropriately used for partial to full-thickness wounds (Stage II - IV), including pressure ulcers, venous ulcers, arterial and diabetic ulcers, donor sites, superficial burn wounds, surgical incisions, abrasions, lacerations and traumatic wounds. This dressing will absorb up to 20 times its own weight in wound exudate. When the dressing comes in contact with wound fluid, it forms a soft absorbent gel-fiber matrix which insulates the wound and creates an ideal moist healing environment. Because of the gel-fiber matrix that is formed, dressing removal is painless and easy. Available in flat dressings or wound packing rope.

## Compares To

- Kaltostat® - Convatec
- DermaGinate – DermaRite
- Algisite – Smith & Nephew
- Maxorb – Medline
- Curasorb® - Kendall
- Sorbsan® - Mylan Bertek

Item No.	Description	Qty.	Measurements
WS7210-A	Alginate Dressing	10/ctn - 10 ctn/cs	2" x 2"
WS7410-A	Alginate Dressing	10/ctn - 10 ctn/cs	4" x 4"
WS71205-A	Alginate Dressing	5/ctn	12" Ribbon



M. C. Johnson Company, Inc.

8801 BUSINESS PARK DRIVE | SUITE 201 | FORT MYERS | FLORIDA | 33912

phone 239.433.0619 | fax 239.433.1498 | www.mcjohnson.com