


BORDERED GAUZE DRESSING

Features and Benefits

- Protects the wound surface
- Soft and flexible
- Gentle on the skin
- Non-adherent permeable contact layer
- Adhesive bordered non-woven cloth tape
- Facilitates wicking of fluid
- Comfortable rounded edges and stretchable fabric tape



Description

This bordered gauze is an absorbent island dressing constructed with a heavy-duty absorbent pad and conformable, non-woven fabric adhesive backing. The porous structure of the Sontara® nonwoven fabric ensures air and water vapor permeability. The fabric is stretchable, conforms well to body contours and allows body movements when in place. The adhesive has been designed to be gentle to the skin but at the same time provide secure fixation. The wound pad can absorb moderate amounts of wound fluids and the low-adherent surface reduces the risk of adherence to the wound.

Compares To

- Telfa Island™ – Kendall
- Medline Bordered Gauze® – Medline
- Covaderm™ – DeRoyal
- Covrsite™ – Smith & Nephew
- Bordered Gauze™ – DermaRite
- Mepore® - MølNycke

Item No.	Description	Qty.	Measurements
WS7425-BGS	Bordered Gauze Dressing	25/ctn	4" x 4" with 2" x 2" Pad
WS7625-BGS	Bordered Gauze Dressing	25/ctn	6" x 6" with 4" x 4" Pad



M. C. Johnson Company, Inc.

8801 BUSINESS PARK DRIVE | SUITE 201 | FORT MYERS | FLORIDA | 33912
phone 239.433.0619 | fax 239.433.1498 | www.mcjohnson.com